

MULTI-FUNCTION AUTOMATIC Blood Pressure Monitor UA-851

Fast. Easy. Accurate.



A multi-function blood pressure monitor in a sleek, contemporary design.

The sleek new Multi-Function Automatic Blood Pressure Monitor combines clinically valuable features with a contemporary new table top design. This monitor includes features such as:

- **Pressure Rating Indicator™** - providing a snapshot of the user's blood pressure classification according to the US JNC7 guidelines.
- **Irregular Heartbeat Detection** - will alert the user to the presence of an irregular heartbeat while still providing an accurate reading.
- **90 memory** – records the last 90 readings.
- Bilingual instruction guide with illustrations and quick start card.
- Four week blood pressure logbook included.

LifeSource blood pressure monitors lead the way in technology and quality – because peoples' lives depend on them. Our commitment to patient-focused thinking has made us the leading primary manufacturer in developing clinically valuable features. The LifeSource blood pressure monitors are fast, easy, accurate and carry a full Lifetime Warranty.

Visit our website at www.LifeSourceOnline.com for more information
or call our Customer Satisfaction hotline at 1-888-726-9966.

- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- 90 Memory Recall
- Displays Average Reading
- Bilingual Instructions
- Four week blood pressure logbook included
- Illuminated start button
- Fast measurement

Multi-Function Automatic Blood Pressure Monitor - UA-851

To obtain an accurate reading, it is important to use the correct cuff size. The UA-851 Multi-Function Automatic Monitor includes either a Medium or a Large cuff. A cuff that is too large will produce a reading that is lower than the correct blood pressure; a cuff that is too small will produce a reading that is higher than the correct blood pressure.

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES	UA-851V	0-93764-60200-9	Medium 9.4-14.2" (24-36cm)
	UA-851VL	0-93764-60201-6	Large 14.2-17.7" (36-45cm)

Each unit comes with a bilingual instruction manual, quick start card, four week blood pressure logbook and warranty card.

Accessories

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES	UA-280	0-93764-04607-0	Medium: 9.4" - 14.2" (24-36cm)
	UA-281	0-93764-01069-9	Large: 14.2" - 17.7" (30 - 45cm)
	TB:181C or TB-233	NA	AC Power Adapter

Technical Specifications

Dimensions	5.9" (w) x 4.9" (l) x 6.1" (h) (150mm x 126mm x 156mm)
Weight	12.3 oz.. (350g) without batteries
Display	Digital, 20-mm character height Pressure/pulse displayed simultaneously
Pressurization	Automatic, using micropump
Deflation	Active exhaust valve
Type	Oscillometric
Accuracy	Pressure: ± 3 mmHg or $\pm 2\%$, whichever is greater Pulse: $\pm 5\%$
Pressure Range	Pressure: 20 mmHg to 280 mmHg Pulse: 40 pulses to 180 pulses
Power Source	4 type "AA" (1.5 volt) alkaline batteries (included) or 120 V AC Adapter (optional)
Battery Life	Approximately 4 months with one daily measurement
Operating Environment	50°F to 104°F (10°C to 40°C) Less than 85% relative humidity
Storage Environment	14°F to 140°F (-10°C to 60°C) Less than 95% relative humidity

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	W x D x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
UA-851V	6.3 x 6.7 x 7.7	2.1	0.19
UA-851VL	6.3 x 7.9 x 7.7	2.1	0.22
Unit Case (10 units)	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
UA-851V	14.2 x 32.3 x 9.3	23.3	2.47
UA-851VL	16.4 x 32.3 x 8.5	23.3	2.61

Ask your A&D Sales Representative about point-of-purchase floor and counter displays.



1756 Automation Parkway
San Jose, CA 95131
www.LifeSourceOnline.com

LIFE SOURCE



Pressure Rating Indicator is a trademark of A&D Medical and is Patent Pending.
LifeSource is a registered trademark of A&D Medical. ©2008 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com