

INSTANT READ Ear Thermometer UT-302

P r e c i s e • I n s t a n t • C o n v e n i e n t



A lightweight and compact thermometer that provides quick and accurate temperature measurements at the push of a button.

The UT-302 Instant Read Ear Thermometer has a unique and attractive design that provides a fast, accurate and easy way to take a temperature through the ear. A temperature reading is provided in just one second! The UT-302 does not require a protective lens cover – its waterproof tip can be easily cleaned with water or rubbing alcohol.

Research indicates that the ear is an ideal site for taking body temperature. The eardrum shares blood vessels with the hypothalamus, the part of the brain that regulates body temperature. Therefore, the ear is an accurate indicator of internal body temperature. Since ear thermometers are designed to receive information from the ear, it is a very safe method for measuring temperature and can be used on children of all ages.

- An accurate reading in just one second (scans the ear 512 times per second)
- Waterproof lens tip that can be easily cleaned with water or rubbing alcohol
- Readings in Fahrenheit or Celsius
- Trilingual instruction guide with illustrations and quick start card in English, French and Spanish

For questions about LifeSource products, call our toll free LifeSource Health Line at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

- Fever Indicator
- 10 Memory Recall
- Illuminated Screen
- Battery Included
- One-button Activation
- One-second Reading
- Trilingual Instructions
- 2 Year Warranty

Instant Read UT-302 Ear Thermometer

- Lightweight unit that provides fast and accurate readings.
- Unique shape allows for a gentle and safe method to obtain a temperature reading.
- Perfect for adults and children of all ages.

PART NUMBER	UPC NUMBER	DESCRIPTION
UT-302	0-93764-60224-5	Ear Thermometer

Each unit comes with a CR-2032 lithium battery, trilingual instruction manual and quick start card. The UT-302 Ear Thermometer is covered by a two-year limited warranty.

Technical Specifications

Model	UT-302
Dimensions	2.2" x 1.75" x 6.2" (55.3 x 44.4 x 158.5 mm)
Weight	2.3 oz. (65.5 g)
Sensor Type	Thermopile Sensor ST118
Measuring range	50°F – 122°F (10°C – 50°C)
Measuring accuracy	±0.4°F (0.2°C): From 96.8°F to 102.2°F (36°C – 39°C) ±0.5°F (0.3°C): Less than 96.8°F (36°C) or greater than 102.2°F (39°C)
Display resolution	0.1°F or °C
Response time	1.0 Second
Memory function	10 memory recall
Thermometer tip	Equipped with a waterproof lens. No lens cover is required.
Automatic shut down	The thermometer will automatically turn off after one minute of non-use.
Power supply	One 3 Volt CR-2032 Lithium battery
Battery life	Approximately 5,000 measurements
Operating environment	60.8°F – 95°F (16°C – 35°C) with relative humidity up to 95%
Storage Environment	-4°F – 122°F (-20°C – 50°C)

Packaging and Shipping Specifications

<u>Selling Unit (Clamshell)</u>	<u>W x D x H (inches)</u>	<u>Wt (lb.)</u>	<u>Individual Item Cube (cubic feet)</u>
UT-302	5.9 x 2.0 x 9.1	0.40	0.06
<u>Unit Case (10 units)</u>	<u>W x D x H (inches)</u>	<u>Wt (lb.)</u>	<u>Case Cube (cubic feet)</u>
UT-302	6.2 x 22 x 9.7	5.0	0.77



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers



1756 Automation Parkway
San Jose, CA 95131
www.LifeSourceOnline.com

LIFE SOURCE®



LifeSource is a trademark of A&D Medical. ©2008 A&D Medical. All rights reserved.

LIT:LS302 11/08

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com